MEDIFLOW® ULTRA II

GCE HEALTHCARE

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GCE[®]
Healthcare

ENGLISH

INSTRUCTION FOR USE: MEDIFLOW® ULTRA II

EN

1. FOREWORD

GCE Medical Regulators are medical devices classified as class IIa according to the Medical Device Directive 93/42/EEC.

Their Compliance with essential requirements of 93/42/EEC Medical Device Directive is based upon EN 10524-4 standard.

This instruction for use covers essential information for the whole life cycle of the regulators:

- Operation
- Cleaning
- Maintenance
- Disposal

2. INTENDED USE

MEDIFLOW® ULTRA is a low pressure regulator with flow-selector device intended for the administration of the following medical gases in the treatment, management, diagnostic evaluation and care of patients:

- oxygen;
- oxygen-enriched air;
- nitrous oxide;
- medical air;
- synthetic air;
- carbon dioxide;
- helium;
- xenon;
- specified mixtures of the gases listed above.

The product is not intended to be used with air or nitrogen for driving surgical tool.

This low pressure regulator is intended to be fitted to medical gas pipeline system terminal units in hospital or ambulance car or to quick connector outlet of medical regulator, with nominal pressures up to 800 kPa.

3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

 **Keep away the product and its associated equipment from:**

- heat sources (fire, cigarettes, etc.),
- flammable materials,
- oil or grease
- water
- dust

 **The product and its associated equipment must be prevented from falling over.**

 **Use only the product and its associated equipment in well vented area.**

 **Always maintain oxygen cleanliness standards.**

Before initial use the product shall be kept in its original packaging. GCE recommends use of the original packaging (including internal sealing bag and caps) if the product is withdrawn from operation (for transport and storage).

Statutory laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-20/+60 °C		-30/+60 °C
	10/100%		10/100%
	600/1200 mbar		600/1200 mbar

4. PERSONNEL INSTRUCTIONS

The Medical Devices Directive 93/42/EEC states that product provider must ensure that all personnel handling the product are provided with the operating instructions & performance data.

⚠ Do not use the product without properly familiarization of the product and its safe operation as defined in this Instruction for use. Ensure user is aware of particular information and knowledge required for the gas in use.

5. PRODUCT DESCRIPTION

Gas from the medical gas pipeline system entering through the inlet connector, then calibrated in the flow control head and deliver to the user through the flow outlet user connection. Different flows can be selected by the flow control knob on the flow head.

A - INLET QUICK CONNECTOR (QC)

The low pressure regulator is fitted to the medical gas pipeline system terminal unit or to medical regulator by a gas specific male probe.

B - FLOW CONTROL HEAD

Different flows can be selected by the flow control knob on the flow head.

Depending on product variant, the flow control head is fitted with a stop function in between maximum flow and zero position.

C - USER OUTLET

The “flow outlet” user connection is the outlet from the flow control head normally with a push on hose fitting (hose nipple) or threaded (for humidifier).

This outlet is for supplying a gas flow (l/min) at atmospheric pressure directly to a patient, e.g. through a cannula or a facemask.

 **The product is not fitted with a pressure relief valve. The user must ensure that an over pressure protection is in place on the medical gas supply pipeline (see requirement of ISO 7396-1) or on medical regulator (see requirement of ISO 10524-1, 3).**

 **The flow outlet shall not be used for driving any medical device.**

Note also that products colour (specially flow control knob) might not follow any gas colour coding.

6. OPERATIONS

6.1. BEFORE USE

VISUAL INSPECTION BEFORE FITTING

- Check if there is visible external damage to the product (including product labels and marking). If it shows signs of external damages, remove it from service and identify its status.
- Visually check if the product is contaminated; if needed, for the regulator, use the cleaning procedure detailed in this section.
- Check if the service or the disposal time has not been exceeded, using GCE or owner's date coding system. If service or disposal time has been exceeded, remove the low pressure regulator from service & suitably identify its status.
- Set flow control knob on ZERO position - until the flow control knob engages.

 **Remove caps from inlet and/or flow outlet. Keep caps in a safe place for reuse during transport or storage.**

 **The low pressure regulator is dedicated only for use with the gas specified on its labelling. Never try to use for another gas.**

FITTING TO MEDICAL GAS PIPELINE SYSTEM TERMINAL UNITS

The low pressure regulator shall be connected only at specified position. Zero flow (“0”) on the flow control knob shall be positioned up (at 12' clock).

LEAKAGE TEST

After connection to medical gas pipeline system terminal unit or to medical regulator, visually check possible leakages:

- low pressure regulator inlet connection to terminal unit or to medical regulator,
- inlet connector to low pressure regulator body,
- user outlet.

 **If any leakage is detected, use the procedure in chapter 6.3 and return the product to GCE for service.**

FUNCTIONAL TEST

- Check that there is a gas flow at each setting (for instance, by listening for the sound of gas flow or checking presence of bubbles in an humidifier).
- Reset flow control knob on ZERO position.

6.2. CONNECTION & USE OF OUTLET

LIST OF RECOGNIZED ACCESSORIES

To be connected to the flow outlet:

- humidifiers, breathing masks or cannulas, gas savers, nebulizers.

 **Before connecting any accessory or medical device to the low pressure regulator, always check that it is fully compatible with the product connection features & the product performances.**

FLOW OUTLET CONNECTION

 **When connecting any accessory to the flow outlet make sure that it is not connected to the patient before operating the product.**

- Ensure the hose/humidifier is compatible with the flow outlet feature.
- Push the hose onto the regulator flow outlet/outlet of humidifier.
- Ensure the hose/humidifier is well engages.

FLOW SETTING

- Set medical gas flow control knob on one of the available flow rates by turning the flow control knob in the requested position.

 **Always ensure that the flow control knob has engaged and not placed between two settings otherwise the flow selector will not deliver correct flow of medical gas.**

 **Do not try to apply an excessive torque on the flow control knob when it stops on the maximum flow position or in zero position.**

 **Only user and doctor is accountable for method of using at specific patient and for assessment of flow value.**

6.3. AFTER USE

- Before disconnecting low pressure regulator from the medical pipe line system terminal unit or medical regulator, reset flow control knob on ZERO position.

- Vent gas pressure from downstream equipment (i.e. via the connected consumers).
- Remove connections from user outlets.
- Refit protection caps. Before refitting the caps ensure they are clean.

7. CLEANING

Remove general contamination with a soft cloth damped in oil free oxygen compatible soap water & rinse with clean water.

Disinfections can be carried out with an alcohol-based solution (spray or wipes).

If other cleaning solutions are used, check that they are not abrasive and that they are compatible the product materials (including labels) and gas.

-  **Do not use cleaning solutions containing ammonia!**
-  **Do not immerse in water or any liquid.**
-  **Do not expose to high temperature (such as autoclave).**

8. MAINTENANCE

8.1. SERVICE AND PRODUCT LIFE TIME

SERIAL NUMBER AND DATE OF PRODUCTION

Form of nine digit serial number stamped on the product is following:

YY MM XXXXX

YY: year of production

MM: month of production

XXXXX: sequence number

Example: serial number 090300521 shows the product produced in March 2009, with sequence number 521.

SERVICE

GCE recommends that a product Periodic inspection is undertaken every year to check proper functionality of the regulator. It should be done by skilled technician.

It should not be assumed that Periodic Inspection recommended by GCE covers every safety procedure or method required by local regulations or statutory requirements, or that abnormal or unusual circumstances may not warrant or suggest further requirements or additional procedures.

LIFE TIME AND WASTE MANAGEMENT

Maximum life time of the product is 10 years.

At the end of the product's life time, the product must be withdrawn from service. The owner of the device shall prevent the reuse of the product and handle the product in compliance with "Directive of European Parliament and Council 2006/12/ES of 5th April 2006 on waste".

In accordance to Article 33 of REACH GCE, s.r.o. as responsible manufacturer shall inform all customers if materials containing 0.1% or more of substances included in the list of Substance of Very High Concern (SVHC). The most commonly used brass alloys used for bodies and other brass components contain 2-3% of lead (Pb), EC no. 231-468-6, CAS no. 7439-92-1. The lead will not be released to the gas or surrounding environment during normal use.

After end of life the product shall be scrapped by an authorized metal recycler to ensure efficient material handling with minimal impact to environment and health. To date we have no information that indicates that other materials containing SVHC of concentrations exceeding 0.1% are included in any GCE product.

8.2. REPAIRS

REPAIRS

Repairs activities cover the replacement of the following damaged or missing components:

- Quick connector
- Piston

The repairs shall be carried out by GCE authorised person only.

Any product sent back to GCE authorised person for maintenance shall be properly packaged. The purpose of the maintenance has to be clearly specified (repair). For product to be repaired a short description of fault and any reference to a claim number might be helpful.

Some repair activities concerning to the replacement of the damaged or missing components can be carried out by the owner of the product. Following parts can be replaced only:

- Caps
- Flow knob and stickers
- Hose nipple (including o-ring)

 **All labels on the equipment must be kept in good, legible condition by the owner and the user during the entire product life time.**

 **All seals and o-rings must be kept in dry, dark and clean environment by the owner and the user during the entire product life time.**

 **Use only original GCE components!**

9. GLOSSARY

EN

	Consult instruction for use		Suitable for Home care use
	Caution		Suitable for Hospital care use
	Keep away from heat and flammable material		Suitable for Emergency care use
	Keep away from oil and grease		Serial number
	Humidity limit		Catalogue number
	Temperature limit		Batch code
	Keep dry		Fragile, handle with care
	Date of manufacture		Manufacturer
	Use by date		Weight of product
	Inlet parameter		Outlet parameter
P_1	Inlet pressure range	P_2	Outlet pressure
P_4	Max outlet pressure (closing pressure)	Q	Outlet flow
	Service or disposal date The serial number indicates the year product has to undergo the overall maintenance activities or has to be disposed off. Refer to the serial number note to determine the Overall maintenance or disposal.		Atmospheric pressure limitation
			Take back equipment for recycling. Do not dispose into unsorted municipal waste.

10. WARRANTY

The Standard Warranty period is two years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the product manufacture shown on the product).

The standard warranty is only valid for products handled according to Instruction for use (IFU) and general industry good practice and standards.

EN

APPENDIX :

Nr 1- Technical and performance data

Nr 2 - Quick coupling feature and connecting/disconnecting procedure

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Czech Republic © GCE, s.r.o.