MEDIMETER®

FLOW-METERING DEVICES; INSTRUCTION FOR USE

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EN

1. FOREWORD

GCE Flowmeters are medical devices classified as class IIa according to the Medical Device Directive 93/42/EEC.

The Compliance with essential requirements of 93/42/EEC Medical Device Directive is based upon ISO 15002 standard.

2. INTENDED USE

Flowmeter is a device intended for the administration of the following medical gases in the treatment and care of patients:

- oxygen O_2
- medical air Air
- nitrous oxide N_2O

The product is not intended to be use with air or nitrogen for driving surgical tool.

This low pressure regulator is intended to be fitted to medical gas pipeline system terminal units in hospital or ambulance car or to quick connector outlet of medical regulator, with nominal pressures up to 5 bar.

3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

 **Keep the product and its associated equipment away from:**

- heat sources (fire, cigarettes, ...),
- flammable materials,
- oil or grease, (especially be careful if hand cream is used)
- water,
- dust.

 **The product and its associated equipment must be prevented from falling over.**

 **Use only the product and its associated equipment in well ventilated area.**

 **Always maintain oxygen cleanliness standards.**

Before initial use the product should be kept in its original packaging. GCE recommends use of the original packaging (including internal sealing bag and caps) if the product is withdraw from operation (for transport, storage). Statutory laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-20/+60 °C		-30/+60 °C
	10/100 %		10/100 %
	600/1200 mbar		600/1200 mbar

⚠ In case of storage at a temperature below -20 °C, do not operate the flow selector until it has been allowed to increase its temperature to a minimum of -20 °C.

4. PERSONNEL INSTRUCTIONS

The Medical Devices Directive 93/42/EEC states that product provider must ensure that all personnel handling the product are provided with the operating instructions & performance data.

⚠ Do not use the product without properly familiarization of the product and its safe operation as defined in this Instruction for use. Ensure user is aware of particular information and knowledge required for the gas in use.


5. PRODUCT DESCRIPTION

Gas from the medical gas pipeline system is entering through the inlet connector (A). Flow is then adjusted to the requested value by turning handwheel (B), which controls the valve for flow regulation. The adjusted flow is visible on the calibrated tube via the float (C). Gas flows through the tube, when constant flow is adjusted, the float is stabilized in a certain position (a certain flow on the scale). Gas is delivered to the user through the flow outlet user connection (D).

A. INLET QUICK CONNECTOR

The flowmeter is fitted to the medical gas pipeline system terminal unit or to medical regulator by a gas specific male probe for specific gas. There are two basic variants:

- 1 – inlet connector is directly mounted on the flowmeter,
- 2 – connection to the rail, including hose and inlet connector.

B. HANDWHEEL

The adjusted flow is visible on the calibrated tube (l/min) via the float. When constant flow is adjusted, the float is stabilized in a certain position on the scale.

C. TUBE WITH THE FLOAT

The adjusted flow is visible on the calibrated tube (l/min) via the float. When constant flow is adjusted, the float is stabilized in a certain position on the scale.

EN

D. USER OUTLET – OUTLET CONNECTION

The “flow outlet” user connection is usually with connection for hose (hose nipple) or threaded (for humidifier).

This outlet is for supplying a gas flow (l/min) atmospheric pressure directly to a patient, e.g. through a cannula or a facemask.

 **The flow outlet shall not be used for driving any medical device.**

Note also that products colour (especially handwheel) might not follow any gas colour coding.

6. OPERATIONS

6.1. BEFORE USE

VISUAL INSPECTION BEFORE FITTING

- Check if there is visible external damage of the flowmeter (including product labels and marking). If it shows signs of external damages, remove it from service and identify its status.
- Visually check if the product is contaminated; if needed, for the flow selector, use the cleaning procedure detailed in this instruction for use.
- Ensure the valve for flow adjustment is completely closed (turning clockwise). When closing don't use an excessive torque (max. recommended closing torque is 0,4 Nm).
- Ensure that connected hose nipple or humidifier to the outlet has correct dimensions.

 **The flowmeter is dedicated only for use with the gas specified on its labeling. Never try to use for another gas.**

FITTING TO MEDICAL GAS PIPELINE SYSTEM TERMINAL UNITS OR TO MEDICAL REGULATOR

- The flowmeter shall be connected only in a defined position – tube up.

FUNCTIONAL TEST

- Check that there is a gas flow when turning the handwheel (for instance, by listening of the sound of gas flow or checking presence of bubbles in a humidifier).
- Close the valve for gas adjustment (turning the handwheel clockwise). When closing don't use an excessive torque.

6.2. CONNECTION & USE OF OUTLET

LIST OF RECOGNIZED ACCESSORIES

To be connected to the flow outlet:

Humidifier, breathing mask or cannula.

EN

-  Before connecting any accessory or medical device to the flowmeter, always check that it is fully compatible with the product connection features & the product performances and there is no negative influence to outlet parameters of the flow selector.

FLOW OUTLET CONNECTION

-  When connecting any accessory to the flow outlet make sure that it is not connected to the patient before operating the product.
- Ensure the hose/humidifier is compatible with the flow outlet feature.
 - Push the hose onto the flowmeter flow outlet/screw onto the humidifier.
 - Ensure the hose/humidifier is well engages.

FLOW SETTING

- Connect the flowmeter to the source of medical gas.
 - Set requested medical gas flow by the handwheel (turning anticlockwise the flow increase, turning clockwise the flow decrease).
-  Do not apply an excessive torque on the flow control knob when it reached the end (stop) position.
-  Only user and doctor is accountable for method of using at specific patient and for assessment of flow value.
-  Fluctuating inlet pressure and ambient temperature can have influence on accuracy of flow measurement.

6.3. AFTER USE

- Before disconnecting flowmeter from the medical pipe line system terminal unit or medical regulator, close the valve for flow regulation (turning the handwheel clockwise). When closing don't use an excessive torque (max. recommended closing torque is 0,4 Nm).
- Remove connected accessories from user outlet of the flowmeter.
- Disconnect flowmeter from the medical gas pipeline system terminal units or medical regulator.

7. CLEANING

Remove general contamination with a soft cloth damped in oil free oxygen compatible soap water & rinse with clean water. Used cleaning solution shall comply with oxygen cleanliness standards for oxygen devices.

Disinfections can be carried out with an alcohol-based solution (spray or wipes).

If other cleaning solutions are used, check that they are not abrasive and that they are compatible the product materials (including labels) and gas.

 Do not use cleaning solutions containing ammonia!

 Do not immerse in water or any liquid.

 Do not expose to high temperature (such as autoclave).

8. LIFETIME, MAINTENANCE AND SERVICE

8.1. PRODUCT LIFE TIME

EN

SERIAL NUMBER AND DATE OF PRODUCTION

Form of nine digit serial number stamped on the product is following:

YY MM XXXXX

YY: year of production

MM: month of production

XXXXX : sequence number

Example: serial number 050300521 shows the regulator produced in March 2005, with sequence number 521.

PRODUCT LIFE TIME AND WASTE MANAGEMENT

Maximum life time of the product is 10 years.

At the end of the product's life time, the product must be withdrawn from service (max.10 years). The owner shall put in place a relevant procedure to ensure the product cannot be used again and follow "Directive 2006/12/EC of the European Parliament and of the Council of 5 April 2006 on waste".

In accordance to Article 33 of REACH GCE, s.r.o. as responsible manufacturer shall inform all customers if materials containing 0.1% or more of substances included in the list of Substance of Very High Concern (SVHC). The most commonly used brass alloys used for bodies and other brass components contain 2-3% of lead (Pb), EC no. 231-468-6, CAS no. 7439-92-1. The lead will not be released to the gas or surrounding environment during normal use. After end of life the product shall be scrapped by an authorized metal recycler to ensure efficient material handling with minimal impact to environment and health.

To date we have no information that indicates that other materials containing SVHC of concentrations exceeding 0.1% are included in any GCE product.

8.2. SERVICE, REPAIRS AND MAINTENANCE

SERVICE

GCE recommends that a product periodic inspection is undertaken at least once a year to do leakage test and check proper functionality of the flow-meter.

LEAKAGE TEST

- After connection to medical gas pipeline system terminal unit or to medical regulator, visually check possible leakages:

- **Leakage test of the valve for flow adjustment:**

1. Close the valve for flow adjustment (turning the handwheel clockwise). When closing don't use excessive torque (max. recommended closing torque is 0,4 Nm).
2. Connect the flow outlet and hose nipple with a hose, the end of the hose immerse into a jar filled with water.
3. Connect the flowmeter to the source of compressed air max. 5 bar.
4. The test: leave the flow selector 1min under pressure. Number of bubbles shall not be higher than 3 bubbles per 1 minut.
5. When test ready, disconnect the flowmeter from the source of compressed air.

- **External leakage test:**

1. Close the valve for flow adjustment (turning the handwheel clockwise). When closing don't use excessive torque (max. recommended closing torque is 0,4 Nm).
2. Close the outlet connection of the flowmeter using a plug, the plug has to have the same thread as the outlet connector of the flowmeter,
3. Connect the flowmeter to the source of compressed air max. 5 bar.
4. Open slowly the valve for flow adjustment (turning the handwheel anti-clockwise).
5. Possible leakage you will detect by application of soap water on product joins.
6. After test disconnect the flowmeter from the source of compressed air, remove the plug from outlet connection and carefully clean the surface of the flowmeter, see chapter 7 (cleaning).

 **Do not use leak detection fluid contain glycols.**

 **If any leakage is detected, use the procedure in chapter 6.3 and return the product to GCE for service.**

 **The product shall not be disassembled or repaired when connected to the source of the pressure.**

Repair activities cover the replacement of the following damaged or missing components:

- Quick connector (probe),
- flow head,
- closing mechanism.

The repairs shall be carried out by a GCE authorised person only. Any product sent back to a GCE authorized person for maintenance shall be properly packaged. The purpose of the maintenance has to be clearly specified (repair, maintenance). For product to be repaired a short description of fault and any reference to a claim number might be helpful. Some repair activities concerning to the replacement of the damaged or missing components can be carried out by a provider / user of the product.

The following parts can be replaced only:

- caps,
- sticker,
- hose nipple (including o-ring),
- handwheel.

 **All labels on the equipment must be kept in good, legible condition by the provider and user during the entire life time.**

 **Use only original GCE components!**

EN

9. GLOSSARY

	Consult Instruction for use		Suitable for Hospital care
	Caution		Suitable for Emergency care
	Keep away from heat and flammable materials		Suitable for Home care
	Keep away from oil and grease		Catalogue number
	Keep dry!		Batch code
	Fragile, handle with care		Serial number
	Temperature limit		Use by date
	Humidity limit		Date of production
	Atmospheric pressure		Manufacturer
	Inlet parameter		Flow reading on the upper edge of the float
	Outlet parameter	P_1	Inlet pressure

10. WARRANTY

The Standard Warranty period is two years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the product manufacture shown on the product).

The standard warranty is only valid for products handled according to Instruction for use (IFU) and general industry good practice and standards.

EN

APPENDIX:

Nr. 1: Technical and performance data

Nr. 2: Quick coupling feature and connecting / disconnecting procedure.

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© GCE, s.r.o.