MEBERIT[®]

IN CASE OF EMERGENCY

Art. 894 "SNAKE"

CE

Art. 895 "SNAKE PLUS"

IT 3 - 8	Nota Informativa
EN 9 - 14	Information Leaflet
DE 15 - 20	Informationsblatt
FR 21 - 26	Note d'Information
ES 27 - 32	Folletto Informativo

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
Figure 1


Figure 2

A	Support structure
B	Restraining belts
C	Transport handles
D	Valve
E	Aluminium aspiration pump (accessory)
F	ABS aspiration pump (accessory)

1. DESCRIPTION AND COMPONENTS

The SNAKE vacuum mattress is a medical device whose primary scope is to immobilise the patient with suspected fractures of the vertebral column and injured or fractured limbs.

The mattress is made in electro-welded spread plastic material and internally contains polystyrene balls; it has an easy-to-use valve with 4 adaptors and is complete with 4 restraining belts and 8 ergonomic handles.

The vacuum mattress has 14 independent chambers, which allow perfect and quick immobilisation of the patient and also allows a uniform distribution of the polystyrene balls.

2. FEATURES

		894	895
Length	mm	2100	2100
Width	mm	1000	1000
Height	mm	60	60
Weight	kg	8	8,5
Maximum capacity	kg	150	150
Special features		14 independent chambers	

The data stated are not binding for the manufacturer; therefore they can undergo variation without forewarning

3. SAFETY

3.1. DIRECTIVE

The vacuum mattress is manufactured according to Regulation 2017/745 on medical devices. Particularly, the vacuum mattress is considered **class I** device.

3.2. DECLARED AND IMPROPER USE

The product is for professional use in sanitary and private environments, with the aim of allowing temporary immobilisation of some parts of the body or for compatible pathologies and eventual transfer onto suitable devices for movement.

The Manufacturer does not assume any liability for use of the product in the conditions stated below, which are considered improper:

	<ul style="list-style-type: none"> ▪ Use by unqualified staff; ▪ Failure to comply with common immobilisation techniques; ▪ Lack of respect of standards/procedures in force regarding safety by the operators; ▪ The use of the device for pathologies that are not envisioned; ▪ Modifications or tampering not authorised by the manufacturer; ▪ Failure to comply with that stated in the information.
-----------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

3.3. INSTRUCTIONS FOR THE OPERATOR

The product is designed for professional use. As well as necessary technical knowledge, the operators must have specific experience in the assistance sector. In particular, they must have acquired practicalness in the use of immobilisation techniques. They must also be in suitable physical conditions i.e. benefit from the force and muscle co-ordination in order to guarantee the maximum level of safety for the patient wherever necessary.

	<ul style="list-style-type: none"> ▪ The vacuum mattresses are not adapt for permanent immobilisation such as a plaster. ▪ Open wounds must be dressed in a sterile manner before hte application of the vacuum mattress, according to the rules regarding aseptic treatment of wounds. ▪ Improper use of the device can cause injury or damage to the patient
-----------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

4. SETTING TO WORK AND USE INSTRUCTIONS

4.1. DELIVERY

The product is delivered to the purchaser inside a cardboard box. On receipt check that the material delivered is in compliance with the order in all of its parts and that it has not undergone damage during transport. Check that the packaging is integral and has not been opened.

	<p><i>For any anomalies, keep the packaging and contact the shipping company and ME.BER within 48 hours. When no specific agreements are made it is intended that the goods are delivered at the purchaser's risk.</i></p>
-------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

4.2. COMMISSIONING

The vacuum mattress can only be commissioned when there is absolute certainty of the efficiency of all of its components and devices: with this aim, ME.BER. qualifide technicians carries out rigorous controls. It is however recommended to carry out a further preventive control according to that indicated in the table.

Test description	Result	
	Posit.	Negat.
Check that all parts of the mattress and any accessories purchased have no visible defects		
Check that the PVC sheets are not torn or cracked		
Check correct functioning of the valve		
Check that the belts allow quick and safe immobilisation.		

If the test result is positive, the device can be used, on the contrary contact ME.BER. technical Assistance Service immediately.

The functional test must be repeated periodically (see §5 "Maintenance") in order to assure perfect efficiency of the device and maximum safety conditions.

4.3. METHOD OF USE

USING THE VACUUM MATTRESS

Perfect application envisions the intervention of at least 4 operators.

1. place the vacuum mattress near to the patient to be immobilised;
2. place the patient inside the mattress, using the means adapt for movement (e.g. scoop stretcher);
3. to put the mattresses in a vacuum condition, the valve must be closed. Introduce the pump's flexible aspiration pipe, with suitable adaptor, directly onto the suction valve, and aspire the air from the mattress (fig. 3);
4. after having reached sufficient rigidity of the mattress, remove the flexible aspiration pipe from the valve;
5. check that the belts are fixed well and are comfortable for the patient;
6. at this point at least 4 (four) operators can lift the mattress, using the relevant handles, and place the immobilised patient on the most suitable aid for transport (fig. 4).

Figure 3

Adjust tightness of the immobilisation belts during aspiration of the mattress.

Movement operations must be carried out with maximum attention and co-ordination

- *Improper use of the device can cause injury and damage to the patient*
- *Creating an excessive vacuum inside the devices or improper use of the immobilisation belts can cause injury and damage to the patient*


Figure 4

RELEASE OF THE LOW PRESSURE

Before storage it is necessary to completely release the low pressure inside the mattress; in order to carry out this operation:

1. turn the valve in an anti-clockwise direction.

4.4. PRESERVATION

The vacuum mattresses are realised with materials that are resistant to environmental conditions envisioned for normal use, they do not require specific attention. However, they must be kept in a closed environment and be protected from dust and dirt in order to guarantee their hygienic conditions. Blows must also be prevented in order not to damage components and/or cause tears, which would make the damaged part unusable. The mattress must be kept in a place that is easily reached by the operators.

Check component integrity and the hygienic conditions of the place of preservation periodically in order to benefit from maximum efficiency and safety at the time of use.

5. MAINTENANCE

5.1. PERIODIC MAINTENANCE TABLE

In order to maintain the product perfectly efficient, periodical maintenance is indispensable.

Whenever a non-conformity to the original product is identified, put the product out of order without delay until the damaged component has been repaired or replaced.

The frequency of routine maintenance is established depending on the use of the product; the table shows the references to be taken into consideration:

Type of intervention	Frequency(interventions/month)		
	≤ 50	51 ÷ 200	> 200
General inspection and correct functioning of the various devices (§4.2)	3 months	2 months	1 month
Cleaning (see successive point)	CLEANING and DISINFECTION must be carried out after every use		

Improper use and maintenance may cause damage to persons and things. Carry out the periodic maintenance according to the indications given in the manual.

5.2. CLEANING

CLEANING and **DISINFECTION** must be carried out **AFTER EVERY USE**

Use a sponge soaked in water and a detergent solution. Rinse and dry.
Disinfect using **EXCLUSIVELY** alcohol-based disinfectants.

*Do not use products containing aggressive agents (bleach, ammonia, etc.) that could cause irreparable damage to the surface of the material.
Do not use abrasive materials such as brushes, steel wool and knife blades to remove any stains*

Failure to carry out cleaning operations can lead to cross infection.

6. ACCESSORIES

Model	Description
888	Aluminium suction pump
896	ABS suction pump
890-001	Case for vacuum mattress

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ME.BER. Srl Unipersonale