MEBER⁺

IN CASE OF EMERGENCY

Art. 706

Art. 708

IT 4 – 5	<i>Nota informativa</i>
EN 6 – 7	<i>Information Leaflet</i>
DE 8 – 9	<i>Informationsblatt</i>
FR 10 – 11	<i>Note d'Information</i>
ES 12 – 13	<i>Nota Informativa</i>

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1. DESCRIPTION

Soft stretcher models 706/708 is a suitable device for transport of patient in case of first aid and medical emergency.

It is entirely realized in nylon, it is light, strong and also fire-proof. Unattacable by molds and bacteriums, it is liquids-proof and easy to clean.

Its strong and double sewin handles (6 for mod.706 and 8 for mod. 708) make easy the handle and a safe transport of the patient.

If requested ME.BER supplies the comfortable and functional shock-proof plastic case for soft stretchers.

2. FEATURES

Classification	Regulation 2017/745 - Class I medical device
----------------	--

		706	708
Length	mm	1950	
Width	mm	750	
N° of handles		6	8
MAX Load	kg	160	
Temperature range	°C	-5 ÷ +50	

The data given are not binding for the manufacturer, therefore they may be modified without notice

3. PROPER AND IMPROPER USE

The product is a professional one and it is intended for an employment in a medical or private situation, to consent the transport by hand, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible

In any case the manufacturer does not assume any responsibility for improper uses of the product under the conditions listed below.

- **Displacement executed by non qualified staff;**
- **Transport of those patients whose pathologies have been defined non compatible by a medical staff;**
- **Non respect of the operators towards specifications in force about safety;**
- **Transport of patients whose weight exceeds 160 Kg;**
- **Modifications or violations which have not been authorized by ME.BER;**
- **Non respect of what has been specified in this information leaflet.**

4. SHIPMENT

The product is delivered to the buyer into a plastic cover. Control, at the moment of receiving, the conformity of the delivered material to the order and also that the material has not been damaged during the transport, by controlling that the packaging is unimpaired and not open.

In case of anomalies, keep the package and contact the carrier within 48 hours from delivery and, if no other specific agreement exists, ME.BER. considers the travelling goods at purchaser's risk

5. SETTING IN USE

The soft stretcher can be set to work only when there is absolute certainty of all its components' efficiency. To be sure of that, qualified technicians of ME.BER. execute strict checks during construction. However it is recommended, before setting to work the soft stretcher, to execute a further preventive check, and an accurate handle washing following the instructions below:

- using a detergent / disinfectant solution clean accurately with a sponge;
- rinse in running water and dry completely.

To clean use sterile rubber latex gloves.

CLEANING and DISINFECTION must be carried out after every each use and before product storage.

6. CONSERVATION

Since the soft stretcher is made of corrosion-proof materials proper for a normal use at environmental condition it does not require special care. However it is necessary to store it in a closed place to protect it from dust and dirt and to assure its hygienic conditions.

Periodically check the integrity of the product and the hygienic conditions of the storage place, in order to assure the maximum efficiency and safety at the moment of utilization.

7. USE

The product has been made for a professional use; it is recommended the employment of qualified operators who have reached great experience in the field of first aid and They also must be in fit physical conditions because they must avail themselves of such strength and muscular coordination to guarantee the maximum level of safety for the patient.

To operate in conditions of safety and efficiency it is necessary the employment of 2 qualified and trained operators on short side (fig. 1). In particular cases, for example. shifings along the stairs, ramps, difficult surfaces or for patient of excessive load, it is necessary the employment of 2 more auxiliary operators on long side (fig. 2). In any case the control and the coordination of the soft stretcher's movements belong to a qualified operator; the others auxiliary operators must contribute to the movements following the instructions of the responsible operator.

- ***If there is suspicion of spinal trauma for the patient, it is recommended an immediate removal of the patient on a proper equipment, such as spinal strretcher, removing, if necessary, the soft stretcher.***
- ***To shift and to place the patient on the soft stretcher follow the medical techniques for the emergency.***
- ***During the transport all the operators must be turned towards the patient.***

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ME.BER. Srl Unipersonale